



°CureWrap

Cooling Therapy can be easy

Your patient is in need of Cooling Therapy.

Induced hypothermia or fever reduction followed by precise temperature control can improve patient outcome.

°CureWrap offers you the means to achieve that goal. It is an integral part of the innovative CritiCool™ System that effectively cools and precisely controls body temperature in any hospital setting.

°CureWrap is a flexible, three-dimensional heat exchange garment. It is designed to provide your patient with optimal surface coverage, maximizing energy transfer.

With °CureWrap, cooling can be promptly initiated and precisely maintained. Easy to handle, hospital staff can prepare the °CureWrap for cooling prior to patient arrival expediting early initiation of treatment.

Using °CureWrap, you can actively control re-warming and achieve a gradual increase in temperature as desired.

Read more about the CritiCool™ System in our associated literature, or contact one of our customer service representatives today.

Simply Set, Wrap & Cool!



Keeping temperature in mind

°CureWrap product principles:

- Non-invasive approach
- Optimal patient surface coverage
- Unique heat removal approach
- Designed for patient care and accessibility
- One size fits all adults
- Body-shaped 3-dimensional garment for intuitive use
- Flexible and durable
- Wide range of cooling therapy applications

°CureWrap design - features & benefits:

- °CureWrap is a one-piece, body shaped garment that facilitates wrapping of your patient.
- °CureWrap maximizes surface coverage to ensure optimal energy transfer.
- Simultaneous water flow through numerous channels provides superior heat exchange capacity.
- Biocompatible, latex-free and anti-static it suits any hospital setting.
- Secured onto your patient using pre-attached hook & loop strips, it allows easy and intuitive access to your patient.
- Permits 72 hours of treatment, °CureWrap is especially designed to facilitate routine patient care and maintenance.



T e c h n i c a l S p e c i f i c a t i o n s

Use with: CritiCool Control Unit
 Set point temperature range: 30°C - 40°C (86°F - 104°F)
 Material: Patient side: layer of skin-compatible non-woven fabric
 Exterior side: brushed loop fabric
 Storage temperature: 10°C - 27°C (50°F - 81°F)
 Working pressure: 0.4 bar
 Maximum pressure: 0.8 bar
 Safety standards and certifications: IEC 60601-1, IEC 60601-2-35, IEC 60601-1-2,
 UL 60601-1, C 22.2

Ordering Information

Product Description	Model	Patient size	Quantity per box
°CureWrap	3500	Adult	8
Connecting tubes*	200-00147	N/A	1

*Supplied with each CritiCool™ System. May be ordered separately.

H e a d q u a r t e r s



Headquarters:
 MTRE Advanced Technologies Ltd.
 4 Ha-Yarden Street,
 Yavne, P.O. Box 102,
 Rehovot 76100, Israel
 Phone: +972-8-9323333
 Fax: +972-8-9328510



United States:
 Mennen Medical Corp.
 101 Witmer Rd.,
 Suite # 100
 Horsham, PA 19044, US
 Phone: (215) 259 1020
 Fax: (215) 675 6212



United Kingdom:
 Charter Kontron Ltd.
 Unit 18 Avant Business Center
 21 Denbigh Road
 Milton Keynes MK1 1DT, UK
 Phone: +44 1908 646070
 Fax: +44 1908 646030



www.mtre.com